



Gesellschaft für
orthopädische Schienen
mbH



ALMOST INVISIBLE



To apply the grid-finger-splint push the Velcro fastener through the slits at the side. Make certain that the soft side is pointing



upwards! Form a loop (see picture) and guide the finger to be protected through this loop into the grid-finger-splint.

**The grid-finger-splints
are available in 10 sizes
and in two colours
transparent or skin
colour.**

Velcro
fasteners
fixate the
splint to the finger.



gfs@fingersplint.de
www.fingersplint.de
+49 4101 377 963

A TOUCH OF SPLINT



www.fingersplint.de

