

# Procedure and Examination Chairs

## ACCESSORY CATALOG



### Models

S2 S3

Sonesta manufactures and sells tables and chairs for urology, gynecology, urodynamic diagnostics and upper GI procedures.



### Leg & foot rests

#### Leg rest holder (single)

This leg rest holder is standard on all models. By loosening one handle, you can adjust the leg rest to any height, width or angle. Can be rotated 360 degrees.

#### Heavy duty rissler leg rests (pair)

This leg rest has a ball and socket construction Wallowing tilt of the leg cushion in any direction. In addition, it has infinite angling possibilities when mounted in our leg rest holder model 6338. Made of stainless steel.

#### Fold up footrests (pair)

These footrests are adjustable in length from 30-42 cm/11,8"-16,5". They are also adjustable in width and slide outward 19 cm/7,5" each. The bottom is in the same upholstery as your chair. When folded they function as an extension cushion that transforms your chair into a table where the patient can lie down. The gas piston easily adjusts the angle. They are easily removed with a click-off system.

#### Footplates (single)

Made of stainless steel with a plastic foot section, fits into our standard leg rest holder.





S3





Art. number

525-6338

525-6336



#### Benefits & use:

- Multiple positioning options
- Patient comfort •
- Enables optimal examination access
- Engineered ergonomics



Standard configuraton

Optional accessory



### Seat cushions

#### Seat Cushion with 170 mm cutout

A 'quick change' accessory seat cushion, with a 170 mm/6,7" deep cutout, that can be switched in seconds with any other seat cushion. You simply lift the front edge, remove it and replace with another type of cushion.

#### Seat cushion with 220 mm cutout

A 'quick change' accessory seat cushion, with a 220 mm/8,7" deep cutout, that can be switched in seconds with any other seat cushion. You simply lift the front edge, remove it and replace with another type of cushion.

#### Seat cushion with 80 mm cutout

A 'quick change' accessory seat cushion for pedriatic care, with a 80 mm/3,1" deep cutout, that can be switched in seconds with any other seat cushion. You simply lift the front edge, remove it and replace with another type of cushion.

#### Seat cushion without cutout

Benefits & use:

A 'quick change' accessory seat cushion that can be switched in seconds with any of our other seat cushions. You simply lift the front edge, remove it and replace with another type of cushion.





Standard configuraton

Optional accessory

• Enables optimal examination access

 Multiple seating options Patient comfort

- Easy to use and interchangeble
- Ergonomic and low profile cushions



S3

S2

Art. number

525-6366-3

525-6366-4

525-6366-5

525-6365



#### Benefits & use:

- Multiple positioning options
- Ability to transform your chair into a flat bed
- Patient comfort
- Adjustable in length, 143 cm—193 cm/56"—76"
- User-friendly
- Ergonomic and low profile cushions



Standard configuraton

Optional accessory

SONESTA LINE ACCESSORY CATALOG





#### Benefits & use:

- Practical use
- Functional design
- Easy to use and interchangeble



Standard configuraton

Optional accessory





### Other accessories

S2 S3 Art. number



#### Benefits & use:

- Versatility
- Helpful examination tools
- Easy access



Standard configuraton

Optional accessory

## How to accessorize

### after your and your patient's needs



Arm rests Easy access to the patient from each side.



Leg rests Swivelling leg rests for ultimate patient positioning and comfort.



Foot rests Hydraulic piston helps switch from foot rests to extension cushion.



Foot plates Convenient foot plates for function and comfort.



Neck pillow An ergonomic accessory that enhances the patient comfort.



Puzzle piece An extension cushion

that transforms your chair into a table or flat bed.



Seat cushion

Choose the perfect fit for your patient and examination.

### Motor control

Ergonomic motor controls, enables precise positioning.



Funnel

A functional funnel kit for performing non-UDS urology procedures.



Paper roll holder

Convenient placement for immediate access the paper towel roll.



Drainage receptacle

A stainless steel receptacle. Perfect for your instruments.



Allows the patient to remain in one place during procedures.

#### An ADDvise Group Company



TF 2018-11-13 Version 4

### Compumed Europe srl

Via di Grotte Portella 6/8 00044 - Frascati (RM) Italia Tel. +39 06 2005441 Fax: +39 06 45200048 email:<u>info@compumedeurope.com</u> web: <u>www.compumedeurope.com</u> SONESTA MEDICAL AB TEGELUDDSVÄGEN 76 115 28 STOCKHOLM, SWEDEN +46 8 502 572 80 INFO@SONESTAMEDICAL.SE ORDER@SONESTAMEDICAL.SE



SONESTA MEDICAL INC. 2 DEBUSH AVENUE UNIT C3 MIDDLETON, MA 01949, USA +1 630-519-3450 INFOUS@SONESTAMEDICAL.SE ORDERUS@SONESTAMEDICAL.SE